


MAXITHERM[®]

Heavy Duty Reusable Vinyl Circulator Pads


Instructions for connecting the Maxitherm[®] to your Gaymar[®] and most other circulators.

Step 1: Remove the old pad from your circulator hose. Remove the Gaymar[®] connectors (both male & Female) from the hose by immersing the connectors into hot water. It may be necessary to cut the connectors off of the hose (save the connectors for future use).



Step 2: The Maxitherm[®] Pad has 2 male connectors intended for the Microtemp[®] circulator. It will be necessary to remove these connectors. Each Maxitherm[®] contains a small package with 4 extra connectors. The two female parts can be inserted into the Gaymar[®] hose while the two male parts will replace the connectors on the Maxitherm[®] pad.

Sm	12" X 18"
Med	22" X 30"
Lg	24" X 60"



Step 3: At this time, your new Maxitherm[®] circulator pad will connect directly to your Gaymar[®] hose and circulator.

* Gaymar[®] is the registered trademark of GAYMAR INDUSTRIES INC.